## INSTRUCTIONS FOR USE

**Macroplastique® Implants (MPQ-1.5 & MPQ-2.5)**  
Used for the Treatment of Stress Urinary Incontinence  
and Vesicoureteral Reflux

STERILE	R
---------	---

### DESCRIPTION

Macroplastique is a soft tissue-bulking agent consisting of solid, irregularly textured, medical grade polydimethylsiloxane elastomer implants suspended in a hydrogel carrier of polyvinylpyrrolidone (PVP or povidone). Macroplastique Implants are held in place at the implantation site when the hydrogel carrier is replaced by body fluids and host fibroblasts subsequently deposit collagen around the implants. The hydrogel is removed by the reticuloendothelial system and excreted unmetabolized from the body through the kidneys.

Histopathologic studies in animals demonstrate an initial inflammatory and subsequent localized cellular response to the implants resulting in the removal of the hydrogel carrier and its replacement by host connective tissue. The initial reaction is characterized by an influx of multinucleated cells and fibroblasts. The reaction becomes quiescent over time, the number of multinucleated cells decrease, fibroblasts mature into fibrocytes, and collagen formation and maturation takes place. Gross examination of excised material reveals discrete implants surrounded by host tissue augmenting (bulking) the original implantation site.

### INDICATIONS

**Stress Urinary Incontinence:** Macroplastique Implants are indicated for the treatment of stress urinary incontinence resulting from an incompetent sphincteric mechanism. Macroplastique may be used as a primary treatment or as a secondary treatment following previous operative procedures. When implanted between the mid-urethra and bladder neck, the implants act as a bulking material aiding in the restoration of urinary continence.

**Vesicoureteral Reflux:** Macroplastique is indicated for the treatment of vesicoureteral reflux. When implanted at the ureteral orifice, the tissue bulking provides occlusive pressure to reduce or eliminate vesicoureteral reflux.

### CONTRAINDICATIONS

Macroplastique is contraindicated in the following conditions:

- acute urogenital tract inflammation
- gross utero-vaginal prolapse
- untreated detrusor instability
- detrusor hyperreflexia
- neuropathic bladder
- overflow incontinence
- pregnancy
- within one year post-partum
- within 12 weeks of a previous Macroplastique treatment
- within 12 weeks of a previous mid-urethral sling placement

**NOTE: Macroplastique is generally not recommended for placement in fibrotic tissue of patients previously treated with polytetrafluoroethylene or rigid particulate implants for incontinence, or in patients with malignancy. At the surgeon's discretion, its use in these conditions may restore continence to an acceptable level.**

### WARNINGS

- Improper patient selection, improper surgical implantation of the implant and/or undercorrection will result in unsatisfactory performance (see "Technique Considerations").
- Injections into the urethra or the subureteric space are procedures requiring precision and caution and should only be performed by surgeons with experience in urological procedures.
- Pre-operative microbiological cultures of urine must be performed to ensure the absence of urinary tract infection.
- Strict aseptic technique must be employed during implantation.

- Use proper endoscopic needle tunneling technique for transurethral injection procedures. See “Endoscopic Procedure-Female Stress Urinary Incontinence” and “Endoscopic Procedure-For the Treatment of Female Stress Urinary Incontinence in Support of Mid-urethral Sling Placement” sections for instructions.
- If the procedure is followed by reactions in excess of those expected (see Adverse Reactions), diagnostic and therapeutic procedures should be implemented at the discretion of the operating surgeon.
- If there is any visual damage to the syringe and/or the inner or outer package, the integrity of the product cannot be guaranteed and should not be used. Return the entire package to Uroplasty or your distributor for replacement.
- Macroplastique is supplied sterile (irradiation sterilized) and is for single use only.
- Do not resterilize.

## ADVERSE REACTIONS

- General complications associated with anesthesia and methods used in the procedure, as well as those associated with the patient’s degree of general tolerance to any implanted foreign material.
- Dysuria, hematuria and frequency of micturition are to be expected postoperatively. If any of these conditions persist past 48 hours, the patient should be instructed to contact the operating surgeon immediately.
- Post-operative retention may occur which may necessitate intermittent catheterization. If the patient remains unable to void freely, continued intermittent catheterization may be necessary. Indwelling urethral catheterization is not recommended.
- Notify Uroplasty or your distributor of product-related complications occurring with the use of this device.

## HOW SUPPLIED

Pre-filled syringes of Macroplastique are supplied sterile and latex-free for single use in two volumes:

MPQ-1.5 Approximately 1.5 ml

MPQ-2.5 Approximately 2.5 ml

Each syringe has a standard male luer lock tip to accommodate the appropriate endoscopic needle and is sealed within a double foil moisture barrier pouch system.

The following implant system components are used with Macroplastique Implants:

Catalog No.	Description
AD	Administration Device
MIS-025	Macroplastique® Implantation System
MIS-015	Macroplastique® Implantation System

Catalog No.	Description
MFN-718	Flexible Endoscopic Needle - Adult
MFN-520	Flexible Endoscopic Needle - Pediatric*
MRN-018	Rigid Endoscopic Needle - Adult
MRN-020	Rigid Endoscopic Needle - Pediatric*

- \* When using pediatric needles, use the MPQ-1.5 syringe rather than the MPQ-2.5 syringe to prevent pressure build-up.

## TECHNIQUE CONSIDERATIONS

- It is recommended that Macroplastique and the endoscopic irrigation fluid temperature be maintained at 25-30°C (77-86°F) during administration. Macroplastique administration may be difficult at lower temperatures.
- When implanting Macroplastique endoscopically, short-beaked, straight or side viewing cystoscopes or nephroscopes are recommended. The use of a resectoscope is not recommended.
- When implanting Macroplastique with the Macroplastique Implantation System (MIS-5.0), refer to the Instructions for Use found in the MIS-025 package in conjunction with this document.
- When implanting Macroplastique with the Macroplastique Sling Support Kit (MIS-SK), refer to the Instructions for Use found in the MIS-015 package in conjunction with this document.

- For best results, inject product slowly. Recommended time to inject 2.5 ml is 2-3 minutes.
- **Note:** The Administration Device (AD) features a rotating bayonet that enables the operator to keep the handle of the Administration Device in a constant position while the needle can be rotated and positioned to align with the desired injection site.
- **Note:** Instructions for attaching needles to either the MPQ-1.5 or MPQ-2.5 syringes. Each syringe of Macroplastique Implants has a standard male luer lock tip to accommodate proper connection to the needles. Be sure to firmly twist and fasten the metal needle hub onto the luer lock tip of the syringe to achieve proper connection.
- **Note:** Only Uroplasty needles with metal hubs can be used with Macroplastique syringes.
- **Note:** The Adapter (ADAPT) must be used to attach a MPQ-1.5 or MPQ-2.5 syringe to the Administration Device (AD).
- Use proper endoscopic needle tunneling technique for transurethral injection procedures. Begin the endoscopic tunneling technique by withdrawing the scope and needle tip into the urethra 1.5 to 2.0 cm distal from the bladder neck. Locate an implantation site and bring the scope up to a 30-45° angle from the surface of the urethral tissue. Puncture the urethral tissue. Advance the needle 0.5 cm in depth. Angle the scope to 0° (parallel with the urethra) and advance the needle another 0.5 cm. Tissue tunneling will limit product loss. Wait approximately 30 seconds prior to withdrawing the needle. See Figures A-D showing proper tissue tunneling.
- Macroplastique placement should be within the submucosal tissue of either the urethra or the subureteric space.
- Superficial implantation of Macroplastique will lead to mucosal damage and loss of the implants. Use the proper endoscopic needle tunneling technique to limit product loss.
- If Macroplastique is injected below the submucosa, the desired tissue bulking will not be achieved. Use the proper endoscopic needle tunneling technique to avoid injecting product too deep.
- Product dries out very quickly. Do not open the package until immediately prior to use. Do not allow needles filled with implants to be exposed to air for an extended period of time. The water portion of the product will evaporate, implants will occlude the needle, and injection will either be difficult or impossible.
- Do not use post-operative indwelling catheters.

### PRE-OPERATIVE CARE FOR ENDOSCOPIC PROCEDURE

Collect and microbiologically culture urine to ensure absence of urinary tract infection. Systemic infection should be considered and ruled out prior to Macroplastique treatment.


1. Prepare the patient with local, regional or general anesthesia.
2. It is recommended to administer systemic antibiotics prior to Macroplastique implantation consistent with current surgical implant procedures. Be certain antibiotics do not interfere with pre-operative microbiological cultures.
3. Place the patient in a modified lithotomy position and surgically prepare per surgeon's instructions.
4. For adults, children and infants, fill bladder to 50% of its capacity.

### ENDOSCOPIC PROCEDURE - FEMALE STRESS URINARY INCONTINENCE (SUI)

1. Peel open the double foil pouch system containing Macroplastique. Remove the syringe from the inner package. Keep the cap on the tip of the syringe. Slide the adapter (ADAPT) over the syringe until it fits down around the syringe flanges. Twist lock the assembled adapter and syringe onto the Administration Device rotating bayonet. Note: Firmly grasp both the adapter and the rotating bayonet to secure the syringe to the Administration Device. The Administration Device features a rotating bayonet that enables the operator to keep the handle of the Administration Device in a constant position while the needle can be rotated and positioned to align with the desired injection site.
2. Using aseptic technique, remove the endoscopic needle from its packaging.
3. Without removing the needle's protective sleeve, attach the endoscopic needle by firmly tightening the metal hub of the needle onto the Macroplastique syringe's Luer lock tip. Remove the needle's protective sleeve. Prior to inserting the needle into the endoscope, express (prime) Macroplastique through it until MACROPLASTIQUE IS VISIBLE AT THE TIP.

- With the patient under local, regional or general anesthesia, advance the endoscope into the bladder. The needle is then passed into the endoscope. When using a cystoscope incorporating an Albarran bridge, the endoscopic needle should be inserted through the bridge with the heel of the needle bevel down to avoid damage to the needle tip. Once the needle has passed and is clear of the bridge assembly, the needle should be rotated 180° with the point of the needle adjacent to the urethral tissue to avoid needle skiing. This needle orientation is recommended in all endoscopic procedures.
- Advance the needle to visualize the needle tip. Withdraw scope and needle tip into the urethra 1.5 to 2.0 cm distal from the bladder neck. Locate the six o'clock position and bring the scope up to a 30-45° angle from the surface of the urethral tissue. Puncture the urethral tissue. Advance the needle 0.5 cm in depth. Angle the scope to 0° (parallel with the urethra) and advance the needle another 0.5 cm. Tissue tunneling will limit product loss (i.e., "toothpaste effect") from the implantation site upon removal of the needle. Correct needle placement within the mucosa is checked by injecting a small amount of the implants. If the needle is properly placed, tissue bulking will immediately appear in the urethral mucosa.

**Figures A-D below show proper tissue tunneling**

- Slowly inject approximately 2.5 ml of Macroplastique at the urethral six o'clock position. After injection is complete, wait approximately 30 seconds before withdrawing the needle from the tissue.
- Repeat the above tunneling and injection procedure at the two o'clock and ten o'clock urethral positions. Minimum recommended Macroplastique volumes at each of these two positions is approximately 1.25 ml per injection site. Product volumes injected may be adjusted to achieve urethral coaptation or closure.
- Depending on the patient's history of a previous surgical incontinence procedure (i.e., bladder neck suspension, sling procedures, etc.), the Macroplastique implantation sites and volumes may be adjusted according to the morphology of the bladder neck and urethra.
- Withdraw the needle prior to removal of the scope from the urethra. Passing the endoscope through the implantation site must be avoided.
- Empty the bladder using a small bore urinary drain (i.e., 10-12 Fr. intermittent catheters). Do not use post-operative indwelling catheters.
- Post-operative care and analgesia are determined by the surgeon. Antibiotic prophylaxis is recommended. Upon discharge, the patient must be able to void freely. In cases of urinary retention, intermittent or suprapubic catheterization is advised until normal voiding resumes. Indwelling urethral catheterization is not recommended.
- If a second implantation of additional material is required, reimplantation should not be done prior to twelve (12) weeks. This allows the initial inflammatory response to subside and the host tissue to fully mature. For this reimplantation, inject distal to the initial Macroplastique placement.


#### **ENDOSCOPIC PROCEDURE - FOR THE TREATMENT OF FEMALE STRESS URINARY INCONTINENCE (SUI) IN SUPPORT OF MID-URETHRAL SLING PLACEMENT**

- Peel open the double foil pouch system containing Macroplastique. Remove the syringe from the inner package. Keep the cap on the tip of the syringe. Slide the adapter (ADAPT) over the syringe until it fits down around the syringe flanges. Twist lock the assembled adapter and syringe onto the Administration Device rotating bayonet. Note: Firmly grasp both the adapter and the rotating bayonet to

secure the syringe to the Administration Device. The Administration Device features a rotating bayonet that enables the operator to keep the handle of the Administration Device in a constant position while the needle can be rotated and positioned to align with the desired injection site.

- Using aseptic technique, remove the endoscopic needle from its packaging.
- Without removing the needle's protective sleeve, attach the endoscopic needle by firmly tightening the metal hub of the needle onto the Macroplastique syringe's Luer lock tip. Remove the needle's protective sleeve. Prior to inserting the needle into the endoscope, express (prime) Macroplastique through it until MACROPLASTIQUE IS VISIBLE AT THE TIP.
- With the patient under local, regional or general anesthesia, advance the endoscope into the bladder. The needle is then passed into the endoscope. When using a cystoscope incorporating an Albarran bridge, the endoscopic needle should be inserted through the bridge with the heel of the needle bevel down to avoid damage to the needle tip. Once the needle has passed and is clear of the bridge assembly, the needle should be rotated 180° with the point of the needle adjacent to the urethral tissue to avoid needle skiing. This needle orientation is recommended in all endoscopic procedures.
- Advance the needle to visualize the needle tip. Withdraw scope and needle tip into the urethra 1.5 to 2.0 cm distal from the bladder neck. Locate the six o'clock position and bring the scope up to a 30-45° angle from the surface of the urethral tissue. Puncture the urethral tissue. Advance the needle 0.5 cm in depth. Angle the scope to 0° (parallel with the urethra) and advance the needle another 0.5 cm. Tissue tunneling will limit product loss (i.e., "toothpaste effect") from the implantation site upon removal of the needle. Correct needle placement within the mucosa is checked by injecting a small amount of the implants. If the needle is properly placed, tissue bulking will immediately appear in the urethral mucosa.

Figures A-D below show proper tissue tunneling

- Slowly inject approximately 1.5 ml of Macroplastique at the 10 o'clock urethral position. After injection is complete, wait approximately 30 seconds before withdrawing the needle from the tissue.
- Repeat the above tunneling and injection procedure at the two o'clock urethral position.
- Depending on the patient's history of a previous surgical incontinence procedure (i.e., bladder neck suspension, mid-urethral sling procedures, etc.), the Macroplastique implantation sites and volumes may be adjusted according to the morphology of the bladder neck and urethra.
- Withdraw the needle prior to removal of the scope from the urethra. Passing the endoscope through the implantation sites must be avoided.
- Empty the bladder using a small bore urinary drain (i.e., 10-12 Fr. intermittent catheters). Do not use post-operative indwelling catheters.
- Post-operative care and analgesia are determined by the surgeon. Antibiotic prophylaxis is recommended. Upon discharge, the patient must be able to void freely. In cases of urinary retention, intermittent or suprapubic catheterization is advised until normal voiding resumes. Indwelling urethral catheterization is not recommended.
- If a second implantation of additional material is required, reimplantation should not be done prior to twelve (12) weeks. This allows the initial inflammatory response to subside and the host tissue to fully mature. For this reimplantation, inject distal to the initial Macroplastique placement.

## **NON-ENDOSCOPIC PROCEDURE - FEMALE STRESS URINARY INCONTINENCE (SUI)**

If Macroplastique is being administered with the Macroplastique Implantation System (MIS-5.0), refer to the Instructions for Use found in the MIS-025 package in conjunction with this document.

## **NON-ENDOSCOPIC PROCEDURE – FOR THE TREATMENT OF FEMALE STRESS URINARY INCONTINENCE (SUI) IN SUPPORT OF MID-URETHRAL SLING PLACEMENT**


If Macroplastique is being administered with the Macroplastique Sling Support Kit (MIS-SK), refer to the Instructions for Use found in the MIS-015 package in conjunction with this document.

## **ENDOSCOPIC PROCEDURE - MALE INCONTINENCE**


The procedure to be used for male incontinence is the same as the endoscopic transurethral procedure for female stress urinary incontinence, with the exception of the following points:

1. Identify the soft, compliant, non-fibrotic tissue in the area of the external sphincter (by reducing irrigation or by requesting the patient to contract the external sphincter).

### **1.a. Functional length of the external sphincter < 1.5 cm.**

### **1.b. Functional length of the external sphincter > 1.5 cm.**

Follow the tissue tunneling injection technique as described in point 6 in the section titled “ENDOSCOPIC PROCEDURE - FEMALE STRESS URINARY INCONTINENCE (SUI)” and in Figures A-D showing proper tissue tunneling.

2. The Macroplastique injection sites should be adjusted according to the morphology of the external sphincter and urethra. Strive for a clockwise positioning with “4 hours” difference between three implantation sites. As a guideline, it is suggested to slowly inject approximately 2.5 ml of Macroplastique at the 6 o’clock position. Wait at least 30 seconds before withdrawing the needle. Slowly inject approximately 1.25 ml to 2.5 ml of Macroplastique at both the 2 o’clock and the 10 o’clock positions. Wait at least 30 seconds before withdrawing the needle. If the patient requires a subsequent treatment, wait at least twelve (12) weeks. This allows the initial inflammatory response to subside, and the host tissue to fully mature.
3. At the time of the subsequent treatment, based on the presence or absence of the prior implant, the Macroplastique injection sites should be adjusted according to the morphology of the external sphincter and urethra. Strive for a clockwise positioning with “4 hours” difference between three implantation sites. At the surgeon’s discretion, slowly implant approximately 1.25 ml to 2.5 ml of Macroplastique 0.5 cm distal of the initial implantation site at each of the 6, 2 and 10 o’clock positions provided the previous implant is no longer present. Wait at least 30 seconds before withdrawing the needle from each site.
4. If further treatment is required, wait for an additional twelve (12) weeks. At the surgeon’s discretion, slowly inject approximately 1.25 ml to 2.5 ml at each of the 2, 10 and 6 o’clock positions, depending on present morphology.
5. Male post-operative care is similar to female post-operative care.

## **ENDOSCOPIC PROCEDURE - VESICoureTERAL REFLUX (VUR)**

**NOTE: When using pediatric needles, use the MPQ-1.5 syringe rather than the MPQ-2.5 syringe to prevent pressure build-up.**

1. Follow steps 1-5 under “Endoscopic Procedure - Female Stress Urinary Incontinence (SUI).”

2. The endoscopic needle should enter the mucosa 0.5 cm anterior to the ureteric orifice. The needle is advanced approximately 0.5 to 1.0 cm staying within the submucosal layer until the needle tip is directly beneath the ureteric orifice. Injecting a small amount of the implants will verify correct needle placement within the mucosa. If the needle is properly placed, tissue bulking will immediately appear at the ureteric orifice.
3. Inject 0.2 to 1.0 ml of Macroplastique into the submucosal layer of the subureteric space until a volcano-like appearance is achieved, resulting in a crescent-shaped ureteric opening located on top of the augmented tissue.
4. Wait approximately 30 seconds before withdrawing the needle from the tissue. This reduces the possibility of product loss (i.e., "toothpaste effect") from the implantation site.
5. Post-operative care is determined by the surgeon. Systemic antibiotics are recommended.
6. If a second implantation of additional material is required, it should not be done prior to twelve (12) weeks. This will allow the initial inflammatory response to subside and the host tissue to fully mature around the implants. For VUR additional treatment, inject beneath the initial Macroplastique placement.

### **PROCEDURE - ADDITIONAL IMPLANTATION (ALL INDICATIONS)**

Prior to reimplantation, it is very important to establish whether or not the initial implantation is still present.

A second implantation may be necessary when insufficient tunneling occurs during the first implantation. When there is insufficient tunneling, the implants may not stay in place. If this is the case, care should be taken to correct the implantation technique. Refer to procedure protocol for specific retreatment instruction for each indication.

### **PATIENT INFORMATION AND INFORMED CONSENT**

Uroplasty relies on the surgeon to advise the patient of all potential risks and complications associated with the proposed surgical procedure, including providing a comparison of the risks, benefits and complications of alternate procedures.

### **PATIENT REGISTRATION FORM**

The Macroplastique Patient Registration Form enclosed within the Macroplastique package must be returned to Uroplasty either directly or through your distributor for product and procedure tracking.

### **DISCLAIMER OF WARRANTIES**

Uroplasty excludes all warranties and responsibilities for:

- improper use of the product,
- failure to follow instructions in this insert (or related Uroplasty inserts, i.e., MIS-025 Instructions for Use or MIS-015 Instructions for Use), and/or,
- failure to return a completed Macroplastique Patient Registration Form.

### **STORAGE CONDITIONS**

Macroplastique should be stored at an ambient temperature of 15-30°C (59-86°F). See product labeling for expiration date.

### **ORDERING INFORMATION**

To order direct or for product information, please contact the Customer Service Department at Uroplasty BV or your local distributor.

Manufacturer:

**Uroplasty BV**

Hofkamp 2

6161 DC Geleen

The Netherlands

Tel.: +31 (0)46 4237920

Fax: +31 (0)46 4237922

email: [customer.service@uroplasty.com](mailto:customer.service@uroplasty.com)

*Macroplastique® is a registered trademark of Uroplasty, Inc. Canadian Patent Number 2133756. EPO Patent Number 0636014. French Patent Number 0636014. German Patent Numbers P3941023.4 & 69318835.9. Italian Patent Number 0636014. Japanese Patent Numbers 2694372 & 3004724. The Netherlands Patent Number 193399. Spanish Patent Number 93908676.5. United Kingdom Patent Numbers 2227176 & 0636014. United States Patent Numbers 5,258,028, 5,571,182 & 5,336,263. Additional Patents Pending.*

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©Uroplasty BV 2007