INSTRUCTIONS FOR USE

PTQ™ Implants (PTQ-2.5)

For the treatment of fecal incontinence due to a weak or disrupted Internal Anal Sphincter.

DESCRIPTION

PTQ™ Implants are solid, irregularly textured, medical grade polydimethylsiloxane elastomer implants suspended in a hydrogel carrier of polyvinylpyrrolidone (PVP or povidone). PTQ Implants are held in place at the implantation site when the hydrogel carrier is replaced by body fluids and host fibroblasts subsequently deposit collagen around the implants. The hydrogel is removed by the reticuloendothelial system and excreted unmetabolised from the body through the kidneys.

Histopathologic studies in animals demonstrate an initial acute inflammatory and subsequent localized cellular response to the implants resulting in the removal of the hydrogel carrier and its replacement by host connective tissue. The initial reaction is characterized by an influx of multinucleated cells and fibroblasts. The reaction becomes quiescent over time, the number of multinucleated cells decrease, fibroblasts mature into fibrocytes, and collagen formation and maturation takes place. Gross examination of excised material reveals discrete implants surrounded by host tissue augmenting (bulking) the original implantation site.

INDICATIONS

Fecal Incontinence: PTQ Implants are indicated for trans-sphincteric implantation in the Internal Anal Sphincter (IAS) for the treatment of fecal incontinence due to a weak or disrupted internal anal sphincter.

CONTRAINDICATIONS

PTQ Implants use in the internal anal sphincter is contraindicated in the following conditions:

- Rectal prolapse
- Fecal impaction
- Symptomatic hemorrhoids
- Perianal and anal scarring
- Perianal sepsis
- Parasitic infestation (thread worms, etc.)
- Congenital anal sphincter defect
- Uncontrolled diabetes
- Immunosuppression
- Acute inflammation, infection, or malignancy
- Pregnancy or within one year postpartum
- Within 12 weeks of a previous PTQ Implants treatment

HOW SUPPLIED

PTQ Implants (PTQ-2.5) are supplied sterile and latex-free for single use in pre-filled syringes containing approximately 2.5 ml of PTQ Implants.

PTQ-2.5 PTQ Implants (2.5 ml)

Each syringe of PTQ Implants has a standard male luer lock to accommodate proper connection to the necessary accessories. The following accessories are recommended for use in the PTQ Implants procedure for soft tissue augmentation:

AD	Administration Device (includes one (1) Adapter)
ADAPT	Adapter for Administration Device
PRN-090	PTQ Rigid Needle

WARNINGS

- PTQ Implants are indicated for trans-sphincteric implantation in the Internal Anal Sphincter (IAS), for the treatment of fecal incontinence due to a weak or disrupted internal anal sphincter. Care should be taken not to implant PTQ Implants intravascularly since potentially damaging foreign body embolisms may result.

- Over-correction of defects must be avoided at all times. A second implantation may be performed at 12 weeks if an additional augmentation is necessary (see “Reimplantation”).
- Improper implantation of the product could result in unsatisfactory performance and may require corrective surgery. The volume of the material implanted is approximately equal to the permanent volume it will occupy when the collagen network is formed (see “Technique Considerations”).
- Strict aseptic technique must be followed during implantation.
- Inspect sealed, sterile package before opening. If either the inner or outer package seal is broken or damaged, the integrity of the product cannot be guaranteed and it should not be used. Return the entire package to your PTQ Implants distributor for replacement.
- Do not open the PTQ Implants package until immediately prior to implantation. Long-term exposure to air will cause dehydration of the product, making it unsuitable for use.
- Do not store the product for reuse once it has been opened. Upon completion of the surgical procedure, discard any unused portion of PTQ Implants.
- Do not resterilize.
- The implantation procedure should not be performed with needles other than those supplied and/or approved by Uroplasty.
- If the procedure is followed by reactions in excess of those that are expected (see “Adverse Reactions”), diagnostic and therapeutic procedures should be implemented at the discretion of the surgeon.
- Pain is usually mild and can be treated with analgesics (see post-operative patient care).

ADVERSE REACTIONS

- Potential complications include all the risks associated with methods utilized in a medical procedure, as well as the patient’s degree of general tolerance to any foreign implanted material.
- Immediately following treatment, transient swelling, bruising, redness and/or discomfort at the implant site may result. Patients should be advised to consult their surgeon if pain, swelling and/or discomfort persists for more than a week.
- Infection may be manifested by fever and/or swelling, tenderness or pain in the surgical area and may appear in either the immediate post-operative period or weeks after implantation. Other reactions may include erosion, undercorrection, lumps and overcorrection, and could necessitate further surgical intervention.

PRE-OPERATIVE PATIENT CARE

- Pre-operative evaluation is essential including complete anorectal physiological testing and endoanal ultrasound.
- Pre-operative laxatives are to be used for 2 days before and 7 days minimum after the procedure (e.g., oral lactulose 15 ml).
- Administer prophylactic broad-spectrum antibiotics immediately prior to implantation (e.g., intravenous injection of 160 mg gentamicin and 500 mg metronidazole, or oral).
- Place the patient in either the Lithotomy or modified “Jack-Knife” (Genu-Pectoral) position. Clean the area with an appropriate surgical solution. Special care should be taken including irrigation of the anal canal with a suitable antimicrobial agent such as Betadine solution.
- Local anesthetic should be used at the site of needle entry and in a track leading up to the area of the internal sphincter to be corrected. Use a local anesthetic with epinephrine to maximize vascular constriction and minimize the risk of hematoma. Infiltrate the skin at the site of needle entry 25 mm from the anal margin with local anesthetic. Continue infiltrating the subcutaneous tissues and in a track leading up to the internal sphincter in the area of the anal cushions. A minimal volume of local anesthetic to achieve optimal anesthesia should be used to avoid tissue distortion.

IMPLANT PROCEDURE - TECHNIQUE CONSIDERATIONS

- To achieve precise, even placement of PTQ Implants, identify the area of sphincteric weakness by performing complete anorectal physiological testing and endoanal ultrasound studies.
- Administer product slowly. If implantation is performed too quickly, uncontrolled amounts of product (PTQ) may form an uncontrolled mass. PTQ Implants can continue to pass through the needle even after the trigger pressure on the Administration Device (Catalog Number AD) is released. After releasing the pressure on the Administration Device, wait approximately 30 seconds before

withdrawal of the needle. Surgeons should familiarize themselves with the flow characteristics of PTQ Implants prior to implantation.

- In the event of excessive bolus formation, massage the area to dissipate the accumulation of product. Avoid over-correction of anal defects.
- Product dries out very quickly. Do not open the package until immediately prior to use. Do not allow needles filled with PTQ Implants to be exposed to air for an extended period of time because the liquid content of the product will evaporate and implants may occlude the needle. If needle blockage occurs, a new needle should be used.
- Store PTQ Implants at an ambient temperature of 15-30° C (59-68° F)
- Do not use expired product. See product labeling for expiration date.

IMPLANT PROCEDURE - DELIVERY SYSTEM PREPARATION

- Peel open the double pouch system, and remove the PTQ Implants syringe from its inner pouch. In order to prevent contamination, DO NOT remove the tip cap before assembly.
- Aseptically remove the PTQ Rigid Needle from its package.
- To prevent damage to the needle tip, DO NOT remove the protective sheath before assembly.
- Slide the adapter (ADAPT) over the syringe until it fits down around the syringe flanges.
- Twist-lock the assembled adapter and syringe onto the Administration Device.
- Remove the protective sheath from the needle. Remove the syringe tip cap and firmly fasten the needle luer lock onto the syringe tip.
- Slide the dosage rod of the Administration Device into the end of the PTQ Implants syringe until the dosage rod contacts the syringe stopper. **NOTE:** The dosage rod may be moved easily, both forward and backward, by depressing the spring-loaded release mechanism and sliding the dosage rod to the desired position.
- Activate the dosage rod on the Administration Device by squeezing the handle until a small amount of PTQ Implants is expressed through the needle. **NOTE:** There is a delay in flow due to the viscosity of PTQ Implants. To stop the flow of product, depress the release mechanism on the Administration Device.

IMPLANT PROCEDURE

- The prepared 18 gauge PTQ Rigid Needle loaded onto the PTQ Implants syringe should enter the skin 25 mm from the anal margin.
- At the level of or slightly above the dentate line, direct the needle medially and upwards, trans-sphincterically (crossing the external sphincter and longitudinal muscle) towards the internal sphincter-submucosal interface. Place a finger in the anal canal to ensure correct direction of injecting needle placement and localization of injected PTQ Implants.
- Care should be taken to avoid penetration of the anal mucosa.
- Slowly inject PTQ Implants. Do not use extra force on the Administration Device dosage rod and gray stopper when the syringe is free of PTQ Implants.
- It is recommended that 2.5 ml of PTQ Implants be injected radially at each of the three positions to correspond with the anatomical locations of the physiological anal cushions in the mid and upper anal canal in the area of the internal anal sphincter (lithotomy approximately 3, 7, and 11 o'clock positions).
- To change to another syringe, disengage the dosage rod by pressing the release mechanism on top of the Administration Device. Retract the dosage rod and remove the syringe and needle. Reload the system as previously directed.
- After the PTQ Implantation procedure is complete, disengage the dosage rod by pressing the release mechanism on top of the Administration Device. Retract the dosage rod and remove the syringe and needle.
- Insertion of instrumentation, such as a proctoscope or ultrasound device, into the anal canal following the procedure is not recommended as bolus disruption and/or loss of tissue bulking may occur.
- Dispose of the implantation needle and used syringes as per standard hospital policy.

POST-OPERATIVE PATIENT CARE

- Patients are discharged after an uneventful post-operative observation with the following medications for 7 days:
 - Oral antibiotics (e.g., oral cephalixin, 500 mg qds and metronidazole, 400 mg tds)
 - Laxatives (e.g., oral lactulose 15 ml)
 - Oral analgesia

ADDITIONAL TREATMENT

In the event of unsatisfactory outcome or continence improvement requiring additional product placement, a retreatment of PTQ Implants should not be carried out until 12 weeks after the initial procedure. Repeating the pre-operative evaluation, including complete anorectal physiological testing and endoanal ultrasound, is essential to identify and investigate the potential reason for unsatisfactory outcome (e.g., incorrect placement) and guide the surgeon with determining the parameters for retreatment.

PATIENT INFORMATION AND INFORMED CONSENT

Uroplasty relies on the surgeon to advise the patient of all potential risks and complications associated with the proposed surgical procedure, including comparison of the risks, benefits and complications of alternate procedures.

PATIENT REGISTRATION FORM

The PTQ Implants Patient Registration Form, which is enclosed in the PTQ Implants product box, must be returned to Uroplasty either directly or through your distributor for product and procedure tracking.

DISCLAIMER OF WARRANTIES

Uroplasty excludes all warranties and responsibilities for:

- Improper use of the product,
- Failure to follow implantation procedures in this insert,
- Failure to return the PTQ Implants Patient Registration Form, and/or,
- Failure to follow the Instructions for Use for any PTQ Implants accessory.

STORAGE CONDITIONS

PTQ Implants should be stored at an ambient temperature of 15-30°C (59-86°F). See product labeling for expiration date.

ORDERING INFORMATION

To order direct or for product information, please contact the Customer Service Department located at Uroplasty BV or your local distributor.

Manufacturer:

Uroplasty BV

Hofkamp 2

6161 DC Geleen, The Netherlands

Tel.: +31 46 423 7920

Fax: +31 46 423 7922

PTQ™ is a trademark of Uroplasty, Inc. U.S. Patent Numbers 5,258,028; 5,336,263 and 5,571,182. United Kingdom Patent Numbers 2,227,176 and 636014. German Patent Numbers DE 3941 023 C2 and DE 693 18 835 T2. The Netherlands Patent Number 193 399. Canadian Patent Number 2133756. Japanese Patent Numbers 2694372 and 3003724. Patents Pending.

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